



Alconox Inc.

Critical Cleaning Experts

30 Glenn St., Suite 309, White Plains, NY 10603 USA

Tel.914.948.4040 • Fax.914.948.4088

www.alconox.com • cleaning@alconox.com

Date COA Created: 6/24/2024

## **CERTIFICATE OF ANALYSIS** **DETERGENT 8**

**Manufactured by: Alconox Inc.**

Sub lot No: 06-17-24-01

Sub lot numbers of Master lot No: E1A2

Manuf. Cat. no.: 1701, 1701-1, 1705, 1715, 1755

Date Of Manufacture: 6/17/2024

Expires end of: 06/26

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

<b><u>ANALYSIS</u></b>	<b><u>RESULT</u></b>	<b><u>ACCEPTED RANGE</u></b>
Appearance:	OK	Clear, water white liquid
pH: (1% solution in DI water)	10.6	10.5 - 11.8
Surface Tension:(1% solution Dyne/cm)	31.1	30.0 - 34.0
Foam (1% in DI water)	NIL	NIL
Specific Gravity:	0.99	0.97 - 1.01

**Laura Falla, QA**

**Alconox, Inc.**

*Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.*