



Alconox Inc.

Critical Cleaning Experts

30 Glenn St., Suite 309, White Plains, NY 10603 USA

Tel.914.948.4040 • Fax.914.948.4088

www.alconox.com • cleaning@alconox.com

Date COA Created: 8/8/2024

## **CERTIFICATE OF ANALYSIS** **DETOJET**

**Manufactured by: Alconox Inc.**

Sub lot No: 07-18-24-01

Sub lot numbers of Master lot No: D1B2

Manuf. Cat. no.: 1601, 1601-1, 1605, 1615, 1630, 1632, 1632-1, 1655

Date Of Manufacture: 7/18/2024

Expires end of: 07/26

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

<b><u>ANALYSIS</u></b>	<b><u>RESULT</u></b>	<b><u>ACCEPTED RANGE</u></b>
Appearance:	OK	Clear, water white liquid
pH: (1% solution in DI water)	11.9	11.5 - 13.0
Foam (1% in DI water)	NIL	NIL
Specific Gravity:	1.34	1.27 - 1.34

**Laura Falla, QA**

**Alconox, Inc.**

*Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.*