



Alconox Inc.

Critical Cleaning Experts

30 Glenn St., Suite 309, White Plains, NY 10603 USA

Tel.914.948.4040 • Fax.914.948.4088

www.alconox.com • cleaning@alconox.com

Date COA Created: 9/17/2024

CERTIFICATE OF ANALYSIS LIQUINOX

Manufactured by: Alconox Inc.

Sub lot No: 09-16-24-01

Sub lot numbers of Master lot No: L1B2

Manuf. Cat. no.: 1201, 1201-1, 1205, 1215, 1230, 1232, 1232-1, 1255

Date Of Manufacture: 9/16/2024

Expires end of: 09/26

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

| <u>ANALYSIS</u> | <u>RESULT</u> | <u>ACCEPTED RANGE</u> |
|--|---------------|-----------------------|
| Appearance: | OK | Clear golden liquid |
| pH: (as is) | 8.2 | 8.0 - 9.5 |
| Surface Tension:(1% solution Dyne/cm) | 28.8 | 27.0 - 31.0 |
| Specific Gravity: | 1.08 | 1.06 - 1.09 |

Laura Falla, QA

Alconox, Inc.

Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.