



Alconox Inc.

Critical Cleaning Experts

30 Glenn St., Suite 309, White Plains, NY 10603 USA

Tel.914.948.4040 ? Fax.914.948.4088

www.alconox.com ? cleaning@alconox.com

Date COA Created: 2/28/2025

CERTIFICATE OF ANALYSIS
DETOJET

Manufactured by: Alconox Inc.

Lot No: 02-25-25-01

IRT Lot No: D1B2

Manuf. Cat. no.: 1601, 1601-1, 1605, 1615, 1630, 1632, 1632-1, 1655

Date Of Manufacture: 2/25/2025

Expires end of: 02/27

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

<u>ANALYSIS</u>	<u>RESULT</u>	<u>ACCEPTED RANGE</u>
Appearance:	OK	Clear, water white liquid
pH: (1% solution in DI water)	12.0	11.5 - 13.0
Foam (1% in DI water)	NIL	NIL
Specific Gravity:	1.34	1.27 - 1.34

Laura Falla, QA

Alconox Inc.

Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.