



Alconox Inc.

Critical Cleaning Experts

30 Glenn St., Suite 309, White Plains, NY 10603 USA

Tel.914.948.4040 • Fax.914.948.4088

www.alconox.com • cleaning@alconox.com

Date COA Created: 4/8/2025

CERTIFICATE OF ANALYSIS **LUMINOX**

Manufactured by: Alconox Inc.

Lot No: 04-04-25-01

IRT Lot No: M1L1

Manuf. Cat. no.: 1901, 1901-1, 1915, 1955

Date Of Manufacture: 4/4/2025

Expires end of: 04/27

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

ANALYSIS

RESULT

ACCEPTED RANGE

Appearance:	OK	Clear to pale yellow liquid
pH: (as is)	6.7	6.0 - 8.0
Surface Tension:(1% solution Dyne/cm)	32.0	32.0 - 38.0
Specific Gravity:	1.03	1.00 - 1.06
Foam (1% in DI water)	NIL	NIL

Laura Falla, QA

Alconox Inc.

***Contains no class 1, class 2, class 3, or class 4 residual
solvent ingredients as defined by ICH Q3C guidelines.***