



Alconox Inc.

Critical Cleaning Experts

30 Glenn St., Suite 309, White Plains, NY 10603 USA

Tel.914.948.4040 • Fax.914.948.4088

www.alconox.com • cleaning@alconox.com

Date COA Created: 5/19/2025

## ***CERTIFICATE OF ANALYSIS*** **ALCOJET**

### **Manufactured by: Alconox Inc.**

Lot No: 05-08-25-11

IRT Lot No: 07-01-24-12

Manuf. Cat. no.: 1401, 1403, 1404, 1404-1, 1425, 1450

Date Of Manufacture: 5/8/2025

Expires end of: 11/28

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

### **ANALYSIS**

### **RESULT**

### **ACCEPTED RANGE**

Appearance:

OK

White powder

Solubility: (1% solution at 25 deg C,  
with agitation)

COMP

Complete

pH:(1% solution in DI water)

11.1

11.0 - 12.5

Surface Tension:(1% solution Dyne/cm)

32.0

31.0-36.0

Foam (1% in DI water)

NIL

NIL

**Laura Falla, QA**

**Alconox Inc.**

***Contains no class 1, class 2, class 3, or class 4 residual  
solvent ingredients as defined by ICH Q3C guidelines.***