



Alconox Inc.

Critical Cleaning Experts

30 Glenn St., Suite 309, White Plains, NY 10603 USA

Tel.914.948.4040 • Fax.914.948.4088

www.alconox.com • cleaning@alconox.com

Date COA Created: 5/27/2025

CERTIFICATE OF ANALYSIS **TERGAJET**

Manufactured by: Alconox Inc.

Lot No: 05-16-25-11

IRT Lot No: T1D3

Manuf. Cat. no.: 2201, 2203, 2204, 2204-1, 2225, 2250

Date Of Manufacture: 5/16/2025

Expires end of: 11/28

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

ANALYSIS

RESULT

ACCEPTED RANGE

Appearance:

OK

White powder

Solubility: (1% solution at 25 deg C,
with agitation)

COMP

Complete

pH:(1% solution in DI water)

11.5

11.0 - 12.5

Surface Tension:(1% solution
Dyne/cm)

32.0

31.0 – 36.0

Foam (1% in DI water)

NIL

NIL

Laura Falla, QA

Alconox Inc.

***Contains no class 1, class 2, class 3, or class 4 residual
solvent ingredients as defined by ICH Q3C guidelines.***