



Alconox Inc.

Critical Cleaning Experts

30 Glenn St., Suite 309, White Plains, NY 10603 USA

Tel.914.948.4040 • Fax.914.948.4088

www.alconox.com • cleaning@alconox.com

Date COA Created: 6/6/2025

CERTIFICATE OF ANALYSIS **DETOJET**

Manufactured by: Alconox Inc.

Lot No: 06-05-25-01

IRT Lot No: 02-26-24-01

Manuf. Cat. no.: 1601, 1601-1, 1605, 1615, 1630, 1632, 1632-1, 1655

Date Of Manufacture: 6/5/2025

Expires end of: 12/28

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

ANALYSIS

RESULT

ACCEPTED RANGE

Appearance:

OK

Clear, water white liquid

pH: (1% solution in DI water)

12.0

11.5 - 13.0

Foam (1% in DI water)

NIL

NIL

Specific Gravity:

1.34

1.27 - 1.34

Laura Falla, QA

Alconox Inc.

Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.