



Alconox, LLC

Critical Cleaning Experts

30 Glenn St., Suite 309, White Plains, NY 10603 USA

Tel.914.948.4040 • Fax.914.948.4088

www.alconox.com • cleaning@alconox.com

Date COA Created: 10/8/2025

## ***CERTIFICATE OF ANALYSIS*** **TERGAJET**

### **Manufactured by: Alconox, LLC**

Lot No: 09-04-25-12

IRT Lot No: T1D3

Manuf. Cat. no.: 2201, 2203, 2204, 2204-1, 2225, 2250

Date Of Manufacture: 9/4/2025

Expires end of: 03/29

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

### **ANALYSIS**

### **RESULT**

### **ACCEPTED RANGE**

Appearance:	OK	White powder
Solubility: (1% solution at 25 deg C, with agitation)	COMP	Complete
pH:(1% solution in DI water)	11.5	11.0 - 12.5
Free alkalinity: (meq/g)	5.1	4.8 - 5.4
Foam (1% in DI water)	NIL	NIL

**Laura Falla, QA**

**Alconox, LLC**

***Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.***