



Alconox, LLC

Critical Cleaning Experts

30 Glenn St., Suite 309, White Plains, NY 10603 USA

Tel.914.948.4040 • Fax.914.948.4088

www.alconox.com • cleaning@alconox.com

Date COA Created: 11/6/2025

## CERTIFICATE OF ANALYSIS ALCOJET

### Manufactured by: Alconox, LLC

Lot No: 10-29-25-11

IRT Lot No: 07-01-24-12

Manuf. Cat. no.: 1401, 1403, 1404, 1404-1, 1425, 1450

Date Of Manufacture: 10/29/2025

Expires end of: 04/29

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

### ANALYSIS

### RESULT

### ACCEPTED RANGE

Appearance:	OK	White powder
Solubility: (1% solution at 25 deg C, with agitation)	COMP	Complete
pH:(1% solution in DI water)	11.4	11.0 - 12.5
Free alkalinity: (meq/g)	4.9	3.7 - 5.5
Foam (1% in DI water)	NIL	NIL

Laura Falla, QA

Alconox, LLC

*Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.*