



Alconox, LLC

Critical Cleaning Experts

30 Glenn St., Suite 309, White Plains, NY 10603 USA

Tel.914.948.4040 • Fax.914.948.4088

www.alconox.com • cleaning@alconox.com

Date COA Created: 1/8/2026

## CERTIFICATE OF ANALYSIS LUMINOX

### Manufactured by: Alconox, LLC

Lot No: 01-06-26-01

IRT Lot No: 04-11-24-01

Manuf. Cat. no.: 1901, 1901-1, 1915, 1955

Date Of Manufacture: 1/6/2026

Expires end of: 07/29

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

### ANALYSIS

### RESULT

### ACCEPTED RANGE

Appearance:	OK	Clear to pale yellow liquid
pH: (as is)	7.1	6.0 - 8.0
APHA color:	0	0 - 150
Specific Gravity:	1.02	1.00 - 1.06
Foam (1% in DI water)	NIL	NIL

Laura Falla, QA

Alconox, LLC

*Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.*