



Alconox, LLC

Critical Cleaning Experts

30 Glenn St., Suite 309, White Plains, NY 10603 USA

Tel.914.948.4040 • Fax.914.948.4088

www.alconox.com • cleaning@alconox.com

Date COA Created: 3/17/2026

**CERTIFICATE OF ANALYSIS**  
**ALCOJET**

**Manufactured by: Alconox, LLC**

Lot No: 03-12-26-11

IRT Lot No: 07-01-24-12

Manuf. Cat. no.: 1401, 1403, 1404, 1404-1, 1425, 1450

Date Of Manufacture: 3/12/2026

Expires end of: 09/29

Storage Conditions: Closed upright and in a cool dry place (should be 5 - 40 deg C or 41 - 104 deg F).

**ANALYSIS**

**RESULT**

**ACCEPTED RANGE**

Appearance:	OK	White powder
Solubility: (1% solution at 25 deg C, with agitation)	COMP	Complete
pH:(1% solution in DI water)	11.4	11.0 - 12.5
Free alkalinity: (meq/g)	4.9	3.7 - 5.5
Foam (1% in DI water)	NIL	NIL

**Laura Falla, QA**

**Alconox, LLC**

*Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.*