



Alconox, LLC

Critical Cleaning Experts

30 Glenn St., Suite 309, White Plains, NY 10603 USA

Tel.914.948.4040 • Fax.914.948.4088

www.alconox.com • cleaning@alconox.com

Date COA Created: 1/6/2026

## ***CERTIFICATE OF ANALYSIS*** **ALCONOX**

### **Manufactured by: Alconox, LLC**

Lot No: 12-17-25-13

IRT Lot No: 03-19-24-14

Manuf. Cat. no.: 1101, 1103, 1104, 1104-1, 1112, 1112-1, 1125, 1150

Date Of Manufacture: 12/17/2025

Expires end of: 06/29

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

### **ANALYSIS**

### **RESULT**

### **ACCEPTED RANGE**

|                                                       |      |                                        |
|-------------------------------------------------------|------|----------------------------------------|
| Appearance:                                           | OK   | White powder with cream colored flakes |
| Solubility: (1% solution at 25 deg C, with agitation) | COMP | Complete                               |
| pH:(1% solution in DI water)                          | 9.2  | 9.0 - 10.0                             |
| Free alkalinity: (meq/g)                              | 1.4  | 1.1 - 2.3                              |

**Laura Falla, QA**

**Alconox, LLC**

***Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.***