



Alconox, LLC

Critical Cleaning Experts

30 Glenn St., Suite 309, White Plains, NY 10603 USA

Tel.914.948.4040 • Fax.914.948.4088

www.alconox.com • cleaning@alconox.com

Date COA Created: 3/26/2026

**CERTIFICATE OF ANALYSIS**  
**TERGAZYME**

**Manufactured by: Alconox, LLC**

Lot No: 03-24-26-12

IRT Lot No: 12-27-23-11

Manuf. Cat. no.: 1301, 1303, 1304, 1304-1, 1325, 1350

Date Of Manufacture: 3/24/2026

Expires end of: 09/29

Storage Conditions: Closed upright and in a cool dry place (should be 5 - 40 deg C or 41 - 104 deg F).

**ANALYSIS**

**RESULT**

**ACCEPTED RANGE**

|   |      |   |
|---|------|---|
| Appearance:   | OK   | White powder with cream colored flakes and small brown granules |
| Solubility: (1% solution at 25 deg C, with agitation) | COMP | Complete  |
| pH:(1% solution in DI water)                          | 9.4  | 9.0 - 10.0  |
| Free alkalinity: (meq/g)                              | 1.1  | 0.8 - 2.6   |

Laura Falla, QA

Alconox, LLC

*Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.*