



Alconox, LLC

Critical Cleaning Experts

30 Glenn St., Suite 309, White Plains, NY 10603 USA

Tel.914.948.4040 • Fax.914.948.4088

www.alconox.com • cleaning@alconox.com

Date COA Created: 4/16/2026

CERTIFICATE OF ANALYSIS
ALCONOX

Manufactured by: Alconox, LLC

Lot No: 04-14-26-11

IRT Lot No: 03-19-24-14

Manuf. Cat. no.: 1101, 1103, 1104, 1104-1, 1112, 1112-1, 1125, 1150

Date Of Manufacture: 4/14/2026

Expires end of: 10 / 29

Storage Conditions: Closed upright and in a cool dry place (should be 5 - 40 deg C or 41 - 104 deg F).

ANALYSIS

RESULT

ACCEPTED RANGE

| | | |
|---|------|--|
| Appearance: | OK | White powder with cream colored flakes |
| Solubility: (1% solution at 25 deg C, with agitation) | COMP | Complete |
| pH:(1% solution in DI water) | 9.1 | 9.0 - 10.0 |
| Free alkalinity: (meq/g) | 1.2 | 1.1 - 2.3 |

Laura Falla, QA

Alconox, LLC

Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.