


ALCONOX, LLC PRODUCT LIST

Available worldwide from leading laboratories, scientific, healthcare and industrial suppliers.

Biodegradable. Concentrated. Interfering residue free.

Alconox, LLC detergents set the standard for critical cleaning. Scientists, engineers and professionals specify Alconox, LLC detergents in their cleaning procedures for optimum economy, environmental safety and performance.

 = spray CIP, machine

 = manual cleaning

 = siphon pipet washer

 = soak, ultrasonic



ALCONOX® Powdered Precision Cleaner

Original general-purpose powder for manual, soak, and ultrasonic cleaning of glass, metal, and plastic. No interfering residues. Has inhibitory residue test.

USDA/FDA compliant. Dilute 1:100. pH 9.5

Size	Mfr. Cat No
4 lb box (1.8 kg)	1104-1
CS 9 x 4 lb (9 x 1.8 kg)	1104
25 lb box (11 kg)	1125
50 lb box (22.5 kg)	1150
100 lb drum (45 kg)	1101
300 lb drum (136 kg)	1103
50 x .5 oz packs (50 x 15g)	1112-1
CS 12 x 50 packs	1112



LIQUINOX® Critical Cleaning Liquid Detergent

Original general-purpose, phosphate-free liquid for manual, soak, and ultrasonic cleaning of glass, metal, and plastic. No interfering residues. Has inhibitory residue test. USDA/FDA compliant. Dilute 1:100. pH 8.5

1 qt bottle (.95 L)	1232-1
CS 12 x 1 qt (12 x .95 L)	1232
1 gal bottle (3.8 L)	1201-1
CS 4 x 1 gal (4 x 3.8 L)	1201
5 gal jerrycan (19 L)	1205
15 gal drum (57 L)	1215
55 gal drum (209 L)	1255



CITRANOX® Liquid Acid Cleaner and Detergent

Liquid acid cleaner for manual, soak and ultrasonic cleaning or acid neutralizing. Corrosion-inhibited for glassware, metals,

and plastic. Surfactant rinse aids for no interfering residues. Has inhibitory residue test. USDA and FDA compliant. Dilute 1:100. pH 2.5

Size	Mfr. Cat No
1 gal bottle (3.8 L)	1801-1
CS 4 x 1 gal (4 x 3.8 L)	1801
5 gal jerrycan (19 L)	1805
15 gal drum (57 L)	1815
55 gal drum (209 L)	1855



DETONOX® Ultimate Precision Cleaner

Non-caustic, liquid concentrate detergent for hand and ultrasonic use on challenging residues. Exceptional removal of cosmetic and pharmaceutical polymers, sticky extracts, resins, creams and lotions. Concentrated, highly emulsifying and penetrating formula is free rinsing for reliable results, and leaves no interfering residues. Dilute 1:100. pH 10.5

1 gal bottle	2301-1
CS 4 x 1 gal (4 x 3.8 L)	2301
5 gal jerrycan (19 L)	2305
55 gal drum (208 L)	2355



LUMINOX® Low-Foaming Neutral Cleaner

Low-foaming, neutral pH detergent for use in manual, ultrasonic, and machine cleaning. Safe-handling detergent used for noncorrosive, light-duty cleaning. No chelation or alkaline waste treat problems. Free rinsing, eliminating interfering residues. Dilute 3:100. pH 7.0

Size	Mfr. Cat No
1 gal bottle (3.8 L)	1901-1
CS 4 x 1 gal (4 x 3.8 L)	1901
15 gal drum (57 L)	1915
55 gal drum (209 L)	1955



TERGAZYME® Enzyme-Active Powdered Detergent

Enzymatic cleaner for manual, soak, and ultrasonic cleaning of proteins, tissue, and body fluids from glass, metals, and plastic, with no interfering residues. Cleans filter membranes, fermenters and bioreactors. Has inhibitory residue test. FDA/USDA compliant. Dilute 1:100. pH 9.5

4 lb box (1.8 kg)	1304-1
CS 9 x 4 lb (9 x 1.8 kg)	1304
25 lb box (11 kg)	1325
100 lb drum (45 kg)	1301
300 lb drum (136 kg)	1303





HALOJET® Low-Foam High Emulsification Cleaner

Low-foaming high emulsifying detergent for superior removal of very tenacious organics, fragrances, dyes found in creams, butters, oils and similar residues. Used in automated, manual, or ultrasonic cleaning. Concentrated and free rinsing CIP formula that does not ship as a hazmat. Dilute 1:100. pH 9.5.

5 gal jerrycan (19L)	2605
55 gal drum (208L)	2655
275 gal tote (1041L)	26275





CITRAJET®

Low-Foaming Liquid Acid Cleaner

Liquid acid cleaner for machine, spray-CIP, manual, and ultrasonic cleaning or acid neutralizing. Corrosion-inhibited

for glassware, metals, and plastic. Surfactant rinse for no interfering residues. Has inhibitory residue test. FDA compliant. Dilute 1:100. pH 2.5

Size	Mfr. Cat No
1 gal bottle (3.8 L)	2001-1
CS 4 x 1 gal (4 x 3.8 L)	2001
5 gal jerrycan (19 L)	2005
55 gal drum (209 L)	2055



DETOJET®

Low-Foaming Liquid Detergent

Heavy-duty alkaline degreaser. Nonionic liquid for labware and machine washers, manual, and

ultrasonic cleaning. Corrosion-inhibited for glass, metals, and plastic. No interfering residues. Has inhibitory residue test. USDA/FDA compliant. Dilute 1:200. pH 12

1 gal bottle (3.8 L)	1601-1
CS 4 x 1 gal (4 x 3.8 L)	1601
5 gal can (19 L) - Air OK	1605
15 gal drum (57 L)	1615
30 gal drum (114 L)	1630
55 gal drum (209 L)	1655



ALCOJET®

Low-Foaming Powdered Detergent

General-purpose nonionic powder for labware and machine washers, manual, and ultrasonic cleaning.

Corrosion-inhibited for glassware, metals, and plastic. No interfering residues. Has inhibitory residue test. USDA authorized. Dilute 1:200. pH 11.5

Size	Mfr. Cat No
4 lb box (1.8 kg)	1404-1
CS 9 x 4 lb (9 x 1.8 kg)	1404
25 lb box (11 kg)	1425
50 lb box (22.5 kg)	1450
100 lb drum (45 kg)	1401
300 lb drum (136 kg)	1403



ALNOCHROMIX™

Oxidizing Acid Additive for Glass Cleaning

Concentrated inorganic persulfate based additive to sulfuric acid to make

a metal-free and surfactant-free glassware cleaner to be used in tissue culture, trace metal, and enzyme kinetics research, etc. Nonmetallic, non-carcinogenic, and easier to dispose replacement for chromic acid.

1 box 10 x 3.1 oz (87 g) pouches	2510
4 boxes of 10 x 3.1 oz (87 g) pouches	2540
3 safety vent caps	2513



TERGAJET®

Low-Foaming Phosphate-Free Powder

Phosphate-free, nonionic powder for machine washers, hand, and ultrasonic cleaning. Corrosion-

inhibited for labware. For water, microbiology, and environmental labs. No interfering residues. Has inhibitory residue test. Dilute 1:200. pH 11.5

Size	Mfr. Cat No
4 lb box (1.8 kg)	2204-1
CS 9 x 4 lb (9 x 1.8 kg)	2204
25 lb box (11 kg)	2225
300 lb drum (136 kg)	2203



SOLUJET®

Low-Foaming Phosphate-Free Liquid

Nonionic liquid for machine, spray-CIP, manual, soak, and ultrasonic cleaning. High emulsifying and penetrating

formula. Corrosion-inhibited for glass, metals, and plastic. No interfering residues. Has inhibitory residue test. FDA compliant. Dilute 1:200. pH 12

1 gal bottle (3.8 L)	2101-1
CS 4 x 1 gal (4 x 3.8 L)	2101
5 gal can (19 L) - Air OK	2105
55 gal drum (209 L)	2155



DETERGENT 8®

Low-Foaming Phosphate-Free Cleaner

Nonionic detergent for PC board washing machines, manual, and ultrasonic cleaning. Ion- and chelate-

free for no conductive residues on electronics and no chelate interference in waste. Phosphate-free. Dilute 3:100. pH 11

Size	Mfr. Cat No
1 gal bottle (3.8 L)	1701-1
CS 4 x 1 gal (4 x 3.8 L)	1701
5 gal jerrycan (19 L) - Air OK	1705
15 gal drum (57 L)	1715



KEYLAJET®

Low-Foaming High Alkaline Liquid

Corrosion inhibited formulation recommended for hard metals, steels, stainless steel, porcelain, ceramic, rigid plastics,

cement, and elastomers. Free rinsing to give you reliable results and no interfering residues. Dilute 1:200. pH 13+

5 gal can (19 L) - Air OK	2405
55 gal drum (208 L)	2455



ALCOTABS®

Critical Cleaning Detergent Tablets

Effervescent detergent tablet for maximum cleaning in siphon pipet washers. Slow release bubbles and

free rinsing give no interfering residues. Has inhibitory residue test. Use 2 tabs/wash. pH 6.5

100 tablet bottle	1500-1
CS 6 x 100 bottles	1500