



Alconox, LLC

Critical Cleaning Experts

30 Glenn St., Suite 309, White Plains, NY 10603 USA

Tel.914.948.4040 • Fax.914.948.4088

www.alconox.com • cleaning@alconox.com

Date COA Created: 5/29/2026

CERTIFICATE OF ANALYSIS
ALCONOX

Manufactured by: Alconox, LLC

Lot No: 05-27-26-11

IRT Lot No: 11-07-25-11

Manuf. Cat. no.: 1101, 1103, 1104, 1104-1, 1112, 1112-1, 1125, 1150

Date Of Manufacture: 5/27/2026

Expires end of: 11 / 29

Storage Conditions: Closed upright and in a cool dry place (should be 5 - 40 deg C or 41 - 104 deg F).

<u>ANALYSIS</u>	<u>RESULT</u>	<u>ACCEPTED RANGE</u>
Appearance:	OK	White powder with cream colored flakes
Solubility: (1% solution at 25 deg C, with agitation)	COMP	Complete
pH:(1% solution in DI water)	9.5	9.0 - 10.0
Free alkalinity: (meq/g)	1.3	1.1 - 2.3

Laura Falla, QA

Alconox, LLC

Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.