



Alconox, LLC

Critical Cleaning Experts

30 Glenn St., Suite 309, White Plains, NY 10603 USA

Tel.914.948.4040 • Fax.914.948.4088

www.alconox.com • cleaning@alconox.com

Date COA Created: 8/17/2026

CERTIFICATE OF ANALYSIS
LEDIZOLV

Manufactured by: Alconox, LLC

Lot No: 08-13-26-01

IRT Lot No: None

Manuf. Cat. no.: 90104

Date Of Manufacture: 8/13/2026

Expires end of: 08/30

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

ANALYSIS

RESULT

ACCEPTED RANGE

Appearance:

OK

pale olive liquid

pH: (as is)

5.5

3.8 - 6.6

APHA color:

30

0 - 200

Laura Falla, QA

Alconox, LLC

Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.