



Alconox, LLC

Critical Cleaning Experts

30 Glenn St., Suite 309, White Plains, NY 10603 USA

Tel.914.948.4040 • Fax.914.948.4088

www.alconox.com • cleaning@alconox.com

Date COA Created: 9/9/2026

CERTIFICATE OF ANALYSIS **DETOJET**

Manufactured by: Alconox, LLC

Lot No: 08-27-26-01

IRT Lot No: 06-05-25-01

Manuf. Cat. no.: 1601, 1601-1, 1605, 1615, 1630, 1632, 1632-1, 1655

Date Of Manufacture: 8/27/2026

Expires end of: 02/30

Storage Conditions: Closed upright and in a cool dry place (should be 15 - 30 deg C or 60 - 90 deg F).

ANALYSIS

RESULT

ACCEPTED RANGE

Appearance:	OK	Clear, water white liquid
pH: (1% solution in DI water)	12.1	11.5 - 13.0
Foam (1% in DI water)	NIL	NIL
Specific Gravity:	1.33	1.31 - 1.36
APHA color:	30	0 - 100

Laura Falla, QA

Alconox, LLC

Contains no class 1, class 2, class 3, or class 4 residual solvent ingredients as defined by ICH Q3C guidelines.