



Alconox, LLC

Critical Cleaning Experts

30 Glenn St., Suite 309, White Plains, NY 10603 USA

Tel.914.948.4040 • Fax.914.948.4088

www.alconox.com • cleaning@alconox.com

Date COA Created: 9/8/2026

## ***CERTIFICATE OF ANALYSIS***

### ***ALCOJET***

#### **Manufactured by: Alconox, LLC**

Lot No: 09-02-26-12

IRT Lot No: 05251421

Manuf. Cat. no.: 1401, 1403, 1404, 1404-1, 1425, 1450

Date Of Manufacture: 9/2/2026

Expires end of: 03 / 30

Storage Conditions: Closed upright and in a cool dry place (should be 5 - 40 deg C or 41 - 104 deg F).

#### **ANALYSIS**

#### **RESULT**

#### **ACCEPTED RANGE**

Appearance:

OK

White powder

Solubility: (1% solution at 25 deg C,  
with agitation)

COMP

Complete

pH:(1% solution in DI water)

11.2

11.0 - 12.5

Free alkalinity: (meq/g)

4.7

3.7 - 5.5

Foam (1% in DI water)

NIL

NIL

**Laura Falla, QA**

**Alconox, LLC**

***Contains no class 1, class 2, class 3, or class 4 residual  
solvent ingredients as defined by ICH Q3C guidelines.***